



Website



YouTube



Facebook



Changzhou Sifary Medcial Technology Co., Ltd

Add: NO.99, Qingyang Road, Xuejia County, Xinbei District, Changzhou City, 213000 Jiangsu, China

Email: ivy@sifary.com

Web: www.eighteeth.com



Landlink GmbH

Add: Dorfstrasse 2/4, 79312, Emmendingen, Germany

Email: info@landlink.eu

DISTRIBUTED BY



CATALOG

E-CONNECT & E-PEX
Cordless Endomotor & Apex Locator

01-06

AirPex
Wireless Apex Locator

07-08

E-connect S
Endomotor With Built-in Apex Locator

09-10

E-xtreme
Make Essence to Extreme

11-12

Ultra X
Cordless Ultrasonic Activator

13-14

Fast-Pack & Fast-Fill
Cordless Obturation System

15-20

UltraMint & UltraMint Pro
All-round Ultrasonic Scaler

21-26

CuringPen
4 LEDs Curing Light

27-28

NanoPix
Intraoral X-Ray Sensor


29-32



E-CONNECT & E-PEX

Cordless Endomotor & Apex Locator





E-CONNECT

Cordless Endo Motor



Torque Calibration



Powerful 1500 mAh Lithium Battery



Digital OLED Display



Integration Option



Low Noise Motor



Screen Cast / Mirroring



Ambidextrous * Display

*User can rotate the Endomotor Screen towards left or right
*In connect function Endomotor mirrors Apex Locator's screen



Five Reciprocation Modes (MO)



Reverse Reciprocation-One Mode
Counter clockwise & clockwise movement

Fwd 30° - Rev 150°
For reciprocating files



Forward Reciprocation-Four Modes
Clockwise & counter clockwise movement

Fwd 150° - Rev 30° Fwd 180° - Rev 30°
Fwd 210° - Rev 30° Fwd 250° - Rev 30°
For all continuous rotation files



Nine Editable Memory Programs (M1 to M9)

RPM Maximum Speed 1000 rpm

Forward Direction of Rotation


360° Clockwise Rotational Movement

Reverse Direction of Rotation

Editable Speed / Rpm Setting

Torque Auto Reverse

Editable Torque setting




One Piece Contra Angle Head

340° Adjustable Handpiece

Detachable Head

Miniature Head

Autoclavable



E-PEX

Smart Apex Locator



Integrate Endomotor with Apex Locator (Connect Function)

- Auto Start & Stop
- Apical Reverse
- Apical Slow Down




Smart Accurate and Reliable Apex Locator

- Advanced Multi - frequency Technology
- High Precision in Wet & Dry Canals
- Large 3.5" LCD Display
- Automatic Calibration
- Compact Size & Stable Design

Screen Display & Its Meaning

	Interface	Interface	Interface	Interface
Display				
Value	12	09	00	-3
Meaning	The corresponding scale value is 2, indicating the distance from the apex is about 2mm	It indicates that the distance from the apex is about 1.5mm	It indicates that the endo file reaches the apex	It indicates that the root canal file has exceeded the apex
Alarm	Didi alarm sound for a long time interval	Didi alarm sound interval becomes shorter	Uninterrupted long beep	Very rapid Didi alarm sound



AirPex

Wireless Apex Locator

International Patent
No.:PCT/CN2019/110429




Feel the Charm of Wireless Charging

Lay it down , Charge it up
Charging is as simple as placing
the AirPex on the charging base



True Wireless Apex Locator

Get rid of the wire to the fullest
Use without magnet wire
Accuracy through wirelessness





E-connect S

Endomotor With Built-in Apex Locator



Torque Calibration
For precise torque output




OLED Screen
Digital OLED Display



Ambidextrous Display
Screen can be rotated as per
Left or right hand use



Powerful Battery
Heavy duty 1500mAh
Lithium battery



M0 - Stand Alone Apex Locator Mode



Apex Locator Bar



Zoom-in & Flash Bar
(when file goes beyond 0.5mm)

APEX LOCATOR

- **Apical control with rotary instruments**
Position of file is displayed on screen and the motor rotates in reverse direction when file reaches pre-set apical working length, to prevent over instrumentation
- **Manual canal measurement**
Manual file clip can be used to measure canals precisely with hand files
- **Multi frequency technology with high accuracy & precision**
- **Apex locator flash bar with zoom-in function**

Working Modes:

FWD - Clockwise Continuous Rotation

REV - Anti-clockwise Continuous Rotation

REC - Reciprocation Movement

- First-in-class cordless open motor with adjustable angles of rotation ranging from 30° to 370° for reciprocation as per user's need

First-in-class feature of speed adjustment in reciprocation mode

ATC - Adaptive Torque Control

ATC manages the file stress in a smart way. When torque limit is reached, instead of going into reverse mode, motor reciprocates till the stress is reduced. Once back to below pre-set torque level, file will again rotate in continuous rotation

Integrated mode features:

- **Auto Start & Stop** - Motor starts working when the file enters inside root canal, and stops when withdrawn
- **Apical Reverse** - When the file reaches pre-set apical reverse point, user can define whether the file starts moving in reverse direction
- **Apical Slow Down** - To preserve the cemento-dentinal junction and avoid file separation, speed of motor slows down when file is in apical region

Specifications:

Rotational Speed Range : 120-1000 rpm
 Torque Range : 0.5-4.0 N.cm
 Memory Programs : 11 (M0 to M10)
 Battery Capacity : 1500 mAh



E-xtreme

Make Essence to Extreme

Smaller | Lighter | Quieter than ever



360° contra-angle rotation



Slim and compact design for maximum comfort and mobility




Miniature contra-angle for better visibility & easier posterior access

– Miniature contra-angle Super Mini head & thin neck



– Normal head & neck




Eighteeth innovative motor, cutting down the noise from inside, allowing you to experience the tranquility.



Less than 100 grams, reducing hand fatigue during treatment, providing the ultimate E-xtreme experience.



Reciprocation function



Built-in file system

Specifications

Model	: E-xtreme	Power supply	: 100 V-240 V
Weight(grams)	: 99 g	Frequency	: 50 Hz / 60 Hz
Torque range	: 0.5 N.cm-4.0 N.cm	Output	: DC 5 V/1 A



Ultra X


Cordless Ultrasonic Activator

Available Tips

Size	Length(mm)
#20 2%	18 & 21
#25 2%	18 & 21



-  Autoclavable Tips
-  Smooth Surface Finish
-  Flexible Tips
*Only X Blue & Silver are Flexible Tips
-  Multiple Tip Options




In Action

Cavitation Effect

It is the process of formation of millions of gas bubbles in the irrigation solution. These gas bubbles grow in size, oscillate in different directions and explode into the liquid, producing powerful cleaning action

Acoustic Streaming

It is a unidirectional flow of liquid produced by ultrasonic waves. Movement of tip produces special effect called acoustic turbulence which accelerates the movement in liquid inside the root canal




Indication

Smear Layer Removal

Effective in removal of smear layer and unplugging of dentinal tubules

Dentin Chips Removal

Ultrasonic agitation forms suspension of residual dentinal chips and tissue remnants with irrigant and facilitates their easy removal



Indication

Biofilm Disruption

Ultrasonic cavitation disrupts the biofilm and helps to reduce bacteria load in the canal

Cleaning of Complex Anatomy

Ultrasonic activation helps irrigants to flush into undetected or un-instrumented canal anatomy like isthmus, fins, lateral canals, apical ramification etc. Also helps to remove vapour lock




Protocol

- Prepare the root canal up to required size with the help of shaping and finishing instruments
- Choose the desired disinfectant and flushing agents for cleaning, debridement and disinfection of canals
- Fill the chosen irrigation solution into the canal, never run Ultra X without irrigating solution
- Select the Ultra X tip of desired size and length which passively fits into canal, always keep 2 - 3mm short of apex
- Activate solution using short vertical up & down strokes for 30 seconds. Be careful to not breach apical 2mm
- Use the device gently and never apply pressure or wedge the tip inside the canal. If any points of resistance, then retrieve the ultrasonic tip
- Replenish the irrigating solution, use suction to remove the loose debris. Repeat the cycle 4 to 5 times
- Completely dry the canal and fill it with choice of obturation material
- It is always recommended to isolate the tooth using rubber dam and use personal protective barriers like head gear, mask, eyeglasses etc








Key Features

-  Ultrasonic Activation Device
-  45 kHz High Frequency
-  Application Based Multiple Tips
-  1500 mAh Powerful Battery







Key Features

-  Two Power Modes with LED Indicator
-  Contra-angle Design
-  Light Weight & Ergonomic
-  Cordless Irrigation Device




Battery

-  1500 mAh Lithium Battery
-  Battery Level Indicator
-  4.5 Hrs of Continuous Operation
-  Option of Direct Plug & Play

Fast-Pack & Fast-Fill

3D Obturation
Devices






Fast-Pack


Obturation Pack Handpiece




Temperature Range

-  - Wide temperature range of 90° to 250°
- Five memory programs to customize temperature and cut-off time as per obturation material used or any other intended use of the device


Internally Heated Tip

-  - Only the working end (4 - 6 mm) of the tip is heated - avoids need for insulating sleeves
- The non-working end is insulated and stays cool avoiding soft tissue burn


Class-leading Controlled Heating

-  - Fast heating - Quick heating prevents PDL damage
- Swift cooling - Swift cooling ensure GP does not follow along the heating tip

Experience Ergonomics

-  - 360° start / stop ring button allows holding the device in any angle comfortable for you
- Left - right handed operation
- Comfortable grip, light in weight



Powerful Battery

-  - Highly efficient 2600mAh Lithium battery
- Direct Plug & Play - In case of low battery, use the device directly with adaptor

Specifications

Heating Time	0.2 seconds for 200 °C 0.5 seconds for 250 °C
Cooling Time	2 seconds for cooling
Charging Time	2 hours to complete charge
Operating Time	4 hours of continuous operation at 250 °C





Protocol for Use

Master Cone Verification & Resizing

Select a master cone which corresponds to the master apical preparation. Check for tug-back and confirm the same with radiograph. Resize the GP by cutting it to 0.5 to 1mm short of the working length. This will prevent peri-apical extrusion of the GP

Selection of Fast-Pack & Plugger Tip

Select the Fast-Pack tip which reaches 4-5 mm short of the working length and give a snug fit at the tip. Mark this length using the silicon stopper. Similarly, select the hand plugger which reaches 4-5mm short of the working length

Gutta Percha & Sealer


Coat the master cone apically with the sealer and insert it into the canal. Brush the canal wall using strokes of GP and sealer

Down Pack with Fast-Pack

Remove the GP from the orifice/coronal part with the help of heated Fast-Pack tip. Plug the GP cone with pre-fit hand plugger at the orifice. Now take the heated Fast-Pack tip 4-5mm short of the working length. The bite of GP will come out of the canal along with the tip. Heated Fast-Pack tip will transfer heat to the GP and convert it into thermoplastic state. Condense the GP with the help of hand plugger to ensure the complete sealing of accessory canals and apex

Fill with Fast-Fill

Filling the canal using the Fast-Fill device in increments of 4mm till the canal is filled completely while condensing the warm GP with the appropriate plugger after each increment




UltraMint & UltraMint Pro

All-round Ultrasonic Scaler

All-round Ultrasonic Scaler


Scaling | Periodontic | Endodontic



Soft Perio Power for Ultimate Comfort

Autoclavable plastic box - tips*6, universal wrench*1, EMS type handpiece*1

Autoclavable plastic box - tips*6, universal wrench*1, Satelec type handpiece*1



*Applicable for all tips from Eighteenth


Smart touch system-ultrasensitive & waterproof



Intelligent self-cleaning-makes disinfection easy and avoids cross-infection




Two ways of water supply- Changing supply methods is easier than ever before



900ml stand-alone bottle with clear 20ml scale mark makes it convenient for formulating liquid.




Detachable handpiece with LED- autoclavable, ergonomic, brighter visibility



Detachable and autoclavable handpiece holder, resists cross-infection



Water control directly from handpiece



Ergonomic foot pedal



• EMS type

• Satelec type

One scaler, 2 solutions

Wide range of frequency: 25kHz-42kHz

– We are proud to give you two solutions to perform at your best—both handpieces (EMS type and Satelec type) provided by Eighteeth are compatible with UltraMint & UltraMint Pro.*
Scaling can be as convenient as you can imagine.

*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.

AI Engine — Scaling has never been so intelligent!



Customized treatment based on tooth surface conditions



Smart touch control system—ultrasensitive & waterproof



Elliptical stroke patterns with polishing effect



Intelligent degenerative feedback




Digital controlled automatic frequency tracking

Ensures machine works on optimum frequency



Sine wave drive

Avoids simulation caused by instantaneous current output; ensures comfortable treatment



• EMS type tips

• Satelec type tips

blue-ringed

• EMS type handpiece


• Satelec type handpiece



Tips: G1 G2 G4 P1 P3 E4 universal wrench, EMS type handpiece

Tips: E3-SD E5-S E7-S P4-S G1-S G2-S universal wrench, Satelec type handpiece

Scaling



Periodontic



Endodontic



Endodontic




Cavity Preparation



*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.

*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.



CuringPen

4 LEDs curing light

Technical parameter:

Dimension: 24*21*230(mm)
Net Weight: 109g
Lens size: 11mm


Power supply:
>1200 times, 3 Sec
Maximum light intensity



CuringPen




Smart controlling chip
keeps light output consistent
Light intensity:
1000mw/cm², 1500mw/cm², 2300mw/cm²
600 mw/cm²(detect mode)



Wavelength:385-515nm
3 Blue, 1 Violet, 10w LED



360° rotatable head




Intensity digital display




Multiple Curing Options

Ramp mode



Pulse mode



Caries detection



CuringPen




NanoPix

Intraoral X-Ray Sensor



1/10 Cost per patient compared to film
3s To complete entire imaging process
50000+ Exposures with no quality loss



Leading technologies : The 1st time in dentistry




AED Technology (Automatic Exposure Detection)

- from chest radiography
- Compatible with almost all DC and AC X-Ray generators
- Digital images comes automatically when your X-Ray is switched on



APS CMOS Technology

- from aerospace industry
- High signal-to-noise ratio with maximum details
- Long life and low power consumption




Direct Deposition CSI Technology

- from chest radiography
- Low dose radiation for safety
- Less exposure time but higher resolution

All-in-one

The most economical X-Ray sensor in the world



Unique design



Ergonomic design

Cupped cable connection ensures comfortable finger positioning



Robust cable connection

70000+ Bendings RA test passed



Ultra-thin body design


4.4mm - the thinnest X-Ray sensor

Easy-of-Use software with only 4 steps


● Step 1. Choose 'Patient'



● Step 2. Image 'Acquisition'






● Step 3. View Image



● Step 4. Generate 'Report'



The Complete Image Solution:
excellent combination of
DR Pad & its base

-  pre-installed software
-  multiple USB ports for diverse needs
-  build an instant connect with users



Specifications(available in 2 sizes)

	NanoPix 1	NanoPix 2
Detector Technology	APS CMOS	APS CMOS
Scintillator	CsI	CsI
Pixel Matrix	1000*1500	1300*1800
Theoretical Resolution	25lp/mm	25lp/mm
True Resolution	20lp/mm	20lp/mm
Active Area (mm)	20*30	26*36
Mechanical Size(mm)	25.4*36.8*4.4	30.4*41.9*4.4
AD Conversion(bit)	16	16
Data Interface	USB 2.0	USB 2.0
Protection	IP68	IP68

Notes:

Notes: